
Thetis Medical


SPLINT SIZING

| | SMALL | LARGE |
|----------|--------------|--------------|
| UK | < 8 | ≥ 8 |
| EU | < 42 | ≥ 42 |
| US men | < 9 | ≥ 9 |
| US women | < 10.5 | ≥ 10.5 |

INSTRUCTIONS

Do not walk in the splint


Do not overtighten the straps


1. Point toes then place shell on leg. Ensure you cover foot but not toes. Fastening patch sits over ankle joint.

2. Lightly do up the leg and foot straps.

3. Position heel in hole of the final strap.

4. Gently fasten.

ENGLISH

BEFORE USING THE SPLINT, PLEASE READ THE FOLLOWING INSTRUCTIONS CAREFULLY TO ENSURE THE SPLINT IS CORRECTLY AND SECURELY POSITIONED.

INTENDED USE: Night-time splint to ensure the foot is kept in plantar flexion (foot pointed down) for the first 3 months after an Achilles tendon tear.

WARNING: Do not stand or try to walk in the splint without crutches. It is not designed to be worn when walking, even on tiptoes. Re-injury to the healing Achilles tendon is possible with misuse of the splint.

TO USE:

- Completely detach the central strap from the splint.
- Place the splint on your lower shin so that it also covers the top of your foot. The foam of the lower end of the splint should sit just above the joints where the toes meet the foot.
- Wrap the upper strap around the back of your lower calf and fasten to itself.
- Wrap the lower strap under your foot and fasten to itself.
- Place your heel in the hole of the central strap. One end of the strap will fasten to the front of the splint and the other end will fasten to the strap itself.

CLEANING: Wipe down surfaces only.

CONTRA INDICATIONS: Use from two weeks post Achilles tear or as instructed by your specialist.

PRECAUTIONS: If you do the straps up too tight you may experience tingling or pain. Loosening the straps will usually resolve this. If there is unusual skin discoloration, swelling or discomfort then please seek medical advice.

For single patient use only.

WHILE EVERY EFFORT HAS BEEN MADE TO ENSURE FUNCTION AND COMFORT, THERE IS NO GUARANTEE THAT RE-INJURY WILL BE PREVENTED WHEN USING THIS SPLINT.

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Thetis Medical Ltd,
Jacon Limited,